Infinity-Lock™ Button System

EN Patient Information Leaflet

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Device Name
Infinity-Lock™ Button System
Product code: 102-1089

Device Type
Orthopaedic Implant – Coracoclavicular Ligament

Device Description
Intended Purpose
The Infinity-Lock Button System is intended to bring the shoulder joint together (Shoulder reduction) whilst it heals following shoulder joint separation.

Intended Performance
The main performance feature of the Infinity-Lock Button System is the ability to reduce the shoulder joint to its normal position and improve the appearance of the shoulder.

Intended Patients
Infinity-Lock Button System is to be used for all patients who are considered suitable for reduction of a separated shoulder joint, with respect to the intended use, indications and contraindications for the device, given in the instructions for use.

Expected Lifetime of the device and necessary follow-up
The Infinity-Lock Button System is a permanent implant intended to remain in the patient’s shoulder for the rest of their life. The device is intended to provide shoulder joint reduction until natural healing can occur, which can be up to approximately two years. The patient should follow the advice of their surgeon who will outline any follow-up requirements after the surgery.

Warnings & Precautions
The Infinity-Lock Button System should not be implanted if it is thought that the patient may have any negative reactions to receiving an implant and appropriate tests should be made before any surgery takes place.

Patients are warned not to exceed appropriate activity levels prescribed by the surgeon or to overload the site of surgery before complete healing has occurred.

The patient and any healthcare professional treating the patient should be aware that Infinity-Lock Button System is compatible with Magnetic Resonance Imaging (MRI scans) with certain conditions. Patients with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, only
- Maximum spatial gradient magnetic field of 4,000 Gauss/cm (40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning (i.e. per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the Infinity-Lock Button System is expected to produce a maximum temperature rise of 1.5° C after 15 minutes of continuous scanning (i.e. per pulse sequence).

In non-clinical testing, the image artifact caused by the Infinity-Lock Button System extends approximately 15 mm from this implant when imaged using a gradient echo pulse sequence and a 3 Tesla MR system.

Information for Safe Use
The patient should follow the rehabilitation programme prescribed by the treating surgeon.

Device Materials
The Infinity-Lock Button System implant contains a tape and a button. The Tape is made from polyethylene terephthalate (Polyester) and colourant D&C Green #6 (less than 0.75 % by weight).

The polyester is composed of the following components:

Component	Concentration (%)
Polyethylene terephthalate	Balance to 100
Fiber finish (water / oil emulsion) (30:70))	≤0.06
Titanium dioxide	0.038-0.042
Diethylene glycol	0.2-0.8
Antimony (Catalyst) traces	<0.02
Aldehyde traces	<0.0001

The button is made from titanium alloy which is composed of the following elements:

Element	Compositional Limits
Mass Fraction %	
Aluminium	5.5 to 6.75
Vanadium	3.5 to 4.5
Iron	0.3 max.
Yttrium	0.005
Oxygen	0.2 max.
Carbon	0.08 max.
Nitrogen	0.05 max.
Hydrogen	0.015 max.
Titanium	Balance to 100

The device also contains two disposable instruments, a Cannulated drill bit and a Guidewire. These instruments are used during the implantation of the device and are made from 100% stainless steel.

Undesirable Side Effects

- Infection
- Pain and stiffness
- Loss of shoulder reduction
- Unstable shoulder joint
- Limited range of motion
- Bone fracture
- Impaired wound healing
- Nerve injuries
- Tissue irritation
- Allergic and/or inflammatory tissue reactions
- Pins and Needles
- Loss of bone

When to contact a health professional
All operations involve some risk. The risk and complications of this procedure are small. Contact a health professional if you experience abnormal levels of the undesirable side effects listed above.

Residual Risks
As with any procedure of this type, there is a risk that surgery may not be effective in treatment or may cause worsening symptoms. New damage or a failure of surgery can occur.

Any serious incident that occurs in relation to the device should be reported to the manufacturer.

A summary of safety and clinical performance can be found at the following link:
<http://ec.europa.eu/tools/eudamed>

Note: The EUDAMED link will only be available after the European database on medical devices, EUDAMED, is launched.

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All rights reserved. © Xiros™ 2024. Worldwide patents and patents pending. Infinity-Lock and Xiros are trademarks of Xiros.
Xiros™ Ltd
Xiros Limited, registered in England No. 1664824.