Poly-Tapes

EN	Patient Information Leaflet
DA	Patientinformationsfolder
DE	Informationsbroschüre für Patientinnen und Patienten
EL	Φύλλο οδηγιών χρήσης για τον ασθενή
ES	Folleto de información para el paciente
ET	Patsiendi infoleht
FI	Pakkausseloste
FR	Notice d'information destinée aux patients
IS	Upplýsingar fyrir sjúklinga
IT	Foglio informativo per il paziente
LT	Informacinis lapelis pacientui
LV	Pacienta informācijas brošūra
NL	Patiënteninformatiefolder
NO	Informasjonsbrosjyre for pasienter
PL	Ulotka informacyjna dla pacjenta
PT	Folheto informativo para o doente
RO	Pliant cu informații destinate pacientului
SL	Navodilo za uporabo
SR	Informativni list za pacijente
SV	Patientinformationsbroschyr
TR	Hasta bilgilendirme broşürü
VI	Tờ Thông tin dành cho Bệnh nhân

Device Name

Poly-Tapes

Product code: 102-1010, 102-1080, 102-1081, 102-1082, 102-1083, 102-1084, 102-1097, 102-1098 & 102-1099

Device Type

Orthopaedic Implant - Tendon/Ligament Reconstruction

Device Description**Intended Purpose**

Poly-Tapes are single use devices intended to be used for the reconstruction of damaged or torn ligaments or tendons.

Intended Performance

The main performance feature of Poly-Tapes is the reconstruction of damaged or torn ligaments or tendons; without the requirement for a tissue graft from the patient, another person or an animal. They help to restore function, reduce pain and allow patients to return to everyday activities, work and sport.

Intended Patients

Poly-Tapes are to be used for all patients who are considered suitable for reconstruction of damaged or torn ligaments or tendons, and that are able to follow standard care instructions following surgery.

Expected Lifetime of the device and necessary follow-up

Poly-Tapes are permanent implants intended to remain in the patient's body for the rest of their life. The device is intended to function while natural healing occurs, which can take up to 5 years. The patient should follow the advice of their surgeon who will outline any follow-up requirements after the surgery.

Warnings & Precautions

Poly-Tapes should not be implanted if it is thought that the patient may have any negative reactions to receiving an implant and appropriate tests should be made before any surgery takes place.

Standard hygiene practices for wound care should be followed.

Patients should not exceed appropriate activity levels prescribed by the surgeon or overload the site of the surgery before healing has occurred.

MRI Safety Information

The patient and any healthcare professional treating the patient should be aware that Poly-Tapes are compatible with Magnetic Resonance Imaging (MRI scans). Patients with Poly-Tapes can be scanned safely in an MR system.

Information for Safe Use

The patient should follow the rehabilitation programme prescribed by the treating surgeon.

Device Materials

Poly-Tapes are made from polyethylene terephthalate (polyester) and meet applicable specifications established by ISO 10993-1 for implantable medical devices.

The polyester is composed of the following components:

Component	Concentration (%)
Polyethylene terephthalate	Balance to 100
Fiber finish (water / oil emulsion) (30:70))	≤0.06
Titanium dioxide	0.038-0.042
Diethylene glycol	0.2-0.8
Antimony (Catalyst) traces	<0.02
Aldehyde traces	<0.0001

The Infinity-Lock tapes contain colorant D&C Green #6 (less than 0.75 % by weight).

Undesirable Side Effects

- Infection
- Pain and stiffness
- Unstable Joint
- Limited range of motion
- Bone fracture
- Impaired wound healing
- Nerve Injuries
- Tissue Irritation
- Allergic and/or inflammatory tissue reactions
- Loss of bone
- Skin breakdown
- Failure of surgery

If any of these side effects occur, Poly-Tapes may need to be removed at the surgeon's discretion.

Poly-Tapes may rupture if prescribed activity levels during rehabilitation are exceeded.

As with any procedure of this type, there is a risk that surgery may not be effective in treatment or may cause worsening symptoms, or new damage can occur.

When to contact a health professional

All operations involve some risk, the risk and complications of this procedure are small. Contact a health professional if you experience abnormal levels of the undesirable side effects listed above.

Residual Risks

As with any procedure of this type, there is a risk that surgery may not be effective in treatment or may cause worsening symptoms.

Any serious incident that occurs in relation to the device should be reported to the manufacturer.

A summary of safety and clinical performance can be found at the following link:
<http://ec.europa.eu/tools/eudamed>

Note: The EUDAMED link will only be available after the European database on medical devices, EUDAMED, is launched.

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Poly-Tapes
Produktkoder: 102-1010, 102-1080, 102-1081, 102-1082, 102-1083, 102-1084, 102-1097, 102-1098 & 102-1099

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All rights reserved. © Xiros™ 2026. Worldwide patents and patents pending. Infinity-Lock and Xiros are trademarks of Xiros.
Xiros™ Ltd
Xiros Limited, registered in England No. 1664824.