Pitch-Patch™ Tissue Reinforcement

EN Patient Information Leaflet

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FI Pakkausseloste

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English

Device Name
Pitch-Patch™ Tissue Reinforcement
Product code: 102-1090XI & 102-1091XI

Device Type
Orthopaedic Implant – Rotator Cuff Reinforcement

Device Description
Intended Purpose
The Pitch-Patch is a single use device intended to be used to reinforce the rotator cuff after or during repair with sutures, if the soft tissue is weak.

Intended Performance
The Pitch-Patch is intended to strengthen the shoulder joint, preserving the rotator cuff repair and the benefits of the repair. The Pitch-Patch is not intended to replace normal body structure or provide the full strength to support the shoulder. Sutures, repairing the injury, and suture anchors, attaching muscle to bone, provide strength for the repair.

The main benefits of Pitch-Patch when used as intended are:

- Strengthening of Rotator-cuff repair
- Lower risk of re-injury
- Improved shoulder function
- Reduced pain
- No risk of foreign tissue rejection, donor disease transmission, poor quality or tissue breakdown as associated with foreign tissue transplants.

The Pitch-Patch has been designed to perform to retain suture strength, which provides stability to the rotator cuff and repair.

Intended Patients
The Pitch-Patch is to be used for patients needing to strengthen the rotator cuff where either the injury cannot be fully repaired using normal methods and / or the soft tissue is weak.

Expected Lifetime of the device and necessary follow-up
The Pitch-Patch is a permanent implant intended to remain in the patient's shoulder for the rest of their life. The device is intended to support a rotator cuff repair where either the injury cannot be completely repaired using normal methods and / or the quality of the soft tissue is weak. The device is intended to provide rotator cuff strength until natural healing can occur, which can be up to two years. The patient should follow the advice of their surgeon who will outline any follow-up requirements after the surgery.

Warnings & Precautions
Patients are warned not to exceed the advised activity levels or to overload the repair before complete healing has occurred.

The Pitch-Patch should not be implanted if it is thought that the patient may have any negative reactions to receiving an implant and appropriate tests should be made before any surgery takes place.

MRI Safety Information
These devices are Magnetic Resonance safe (i.e., an item that poses no known hazards in all MR environments).

Information for Safe Use
The patient should follow the rehabilitation programme prescribed by the treating surgeon.

Device Materials
The Pitch-Patch is made from non-absorbable polyethylene terephthalate (Polyester).

The polyester is composed of the following components:

Component	Concentration (%)
Polyethylene terephthalate	Balance to 100
Fiber finish (water / oil emulsion) (30:70))	≤0.05
Titanium dioxide	0.038-0.042
Diethylene glycol	0.2-0.8
Antimony (Catalyst) traces	<0.02
Aldehyde traces	<0.0001

Undesirable Side Effects

- Inflammation
- Allergic reaction
- Irritation
- Infections
- Wound healing complications
- Frozen shoulder
- Creaking of shoulder joint
- Swelling of the shoulder
- Re-tear

If any of these conditions occur, the device may need to be removed at the surgeon's discretion.

The Pitch-Patch may rupture, or the sutures used to secure the Pitch-Patch may break through the device, or the sutures may break through the soft tissue if prescribed activity levels during rehabilitation are exceeded.

As with any procedure of this type, there is a risk that surgery may not be effective in treatment, may cause worsening symptoms, or new damage may occur.

When to contact a health professional
All operations involve some risk, the risk and complications of this procedure are small. Contact a health professional if you experience abnormal levels of the above listed undesirable side effects.

Residual Risks
As with any procedure of this type, there is a risk that surgery may not be effective in treatment or may cause worsening symptoms. New damage or a failure of surgery can occur.

Any serious incident that occurs in relation to the device should be reported to the manufacturer.

Australia: Any serious incident that occurs in relation to the device should be reported to the manufacturer and the Therapeutic Goods Administration www.tga.gov.au

A summary of safety and clinical performance can be found at the following link: <http://ec.europa.eu/tools/eudamed>

Note: The EUDAMED link will only be available after the European database on medical devices, EUDAMED, is launched.

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All rights reserved. © Xiros™ 2026. Worldwide patents and patents pending. Pitch-Patch and Xiros are trademarks of Xiros.
Xiros™ Ltd
Xiros Limited, registered in England No. 1664824.