JewelACL™

EN	Patient Information Leaflet
DA	Patientinformationsfolder
DE	Informationsbroschüre für Patientinnen und Patienten
EL	Φύλλο οδηγιών χρήσης για τον ασθενή
ES	Folleto de información para el paciente
ET	Patsiendi infoleht
FI	Pakkausseloste
FR	Notice d'information destinée aux patients
IS	Upplýsingar fyrir sjúklinga
IT	Foglio informativo per il paziente
LT	Informacinis lapelis pacientui
LV	Pacienta informācijas brošūra
NL	Patiënteninformatiefolder
NO	Informasjonsbrosjyre for pasienter
PL	Ulotka informacyjna dla pacjenta
PT	Folheto informativo para o doente
RO	Pliant cu informații destinate pacientului
SL	Navodilo za uporabo
SR	Informativni list za pacijente
SV	Patientinformationsbroschyr
TR	Hasta Bilgilendirme Broşürü
VI	Tờ Thông tin dành cho Bệnh nhân

English

Device Name
JewelACL™, 102-6003

Device Type
Orthopaedic Implant – ACL Reconstruction

Device Description
Intended Purpose
The JewelACL is intended to be used for reconstruction of the anterior cruciate ligament.

Intended Performance
The JewelACL is similar in strength to the natural ACL when implanted. It may be implanted on its own, or with a hamstring tendon to reconstruct your ACL.

The expected clinical benefits of the implant are:

- Restoration of knee joint stability
- Reduction in knee joint laxity
- Restoration of full range of motion
- Rapid return to work and sport
- No donor site (when used on its own)
- Lower risk of pain or weakness from donor site (when used with your hamstring tendon)

Intended Patients
JewelACL is for all patients who require ACL reconstruction, with respect to the indications and contraindications. Speak to your clinician for further details.

Expected Lifetime of the device and necessary follow-up
The JewelACL is a permanent implant. It is intended to remain implanted during your whole lifetime. It is not intended to be removed. It is expected to function for up to 2 years. After this time, the natural tissue formed on and around the implant is expected to take over the function.

Warnings & Precautions
There are no data to support the use of the JewelACL in pregnant and breast-feeding women.

You should not exceed the activity levels advised by your clinician.

MRI Safety Information
The JewelACL is MR safe (i.e. an item that poses no known hazards if you have an MRI scan).

Information for Safe Use
You should follow the rehabilitation programme prescribed by the treating clinician.

Device Materials
Polyethylene terephthalate (polyester). The polyester has been subjected to a gas plasma process. This gives the device a high affinity for water so native tissue can integrate into the device as part of a natural body response.

The polyester is composed of the following components.

Component	Concentration (%)
Polyethylene terephthalate	Balance to 100
Fiber finish (water / oil emulsion) (30:70))	≤0.06
Titanium dioxide	0.038-0.042
Diethylene glycol	0.2-0.8
Antimony (Catalyst) traces	<0.02
Aldehyde traces	<0.0001

Potential Undesirable Side Effects

- Infections
- Sensitivity and adverse reaction to implant materials.
- Knee laxity or reduced range of motion

If any of these conditions occur, the device may need to be removed at the surgeon’s discretion.

As with any procedure of this type, there is a risk that surgery may not be effective in treatment or may cause worsening symptoms. New damage or re-rupture may occur.

When to contact a health professional
All operations involve some risk, the risk and complications of this procedure are small. Contact a health professional if you experience abnormal levels of the above listed undesirable side effects.

Residual Risks
As with any procedure of this type, there is a risk that surgery may not be effective in treatment and symptoms may continue or worsen. New damage or re-rupture can occur.

Any serious incident that occurs in relation to the device should be reported to the manufacturer.

Australia: Any serious incident that occurs in relation to the device should be reported to the manufacturer and the Therapeutic Goods Administration www.tga.gov.au

A summary of safety and clinical performance can be found at the following link:
<http://ec.europa.eu/tools/eudamed>

Note: The EUDAMED link will only be available after the European database on medical devices, EUDAMED, is launched.

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All rights reserved. © Xiros™ 2026. Worldwide patents and patents pending. JewelACL and Xiros are trademarks of Xiros.
Xiros™ Ltd
Xiros Limited, registered in England No. 1664824.