




XIROS


Poly-Tapes

For soft tissue approximation
Indicated for patients requiring soft tissue
approximation and reconstruction of
ligaments and tendons

Safe
Strong
Simple

Why Poly-Tapes?

- The open weave structure of **Poly-Tapes** allows a natural body response whereby native tissue integrates into the device
- Made from woven multifilament high tenacity polyester (polyethylene terephthalate) fibre, providing a soft flexible structure with high strength
- Range of sizes available to suit varying patient anatomy
- No requirement to harvest a soft tissue graft to reconstruct the ligament or tendon; reduces morbidity



Order codes

Standard **Poly-Tapes** (supplied sterile):

102-1010 Poly-Tape, 10 mm x 500 mm (open weave)
102-1080 Poly-Tape, 10 mm x 800 mm (open weave)
102-1081 Poly-Tape, 15 mm x 800 mm (open weave)
102-1082 Poly-Tape, 20 mm x 800 mm (open weave)
102-1083 Poly-Tape, 30 mm x 800 mm (open weave)
102-1084 Poly-Tape, 40 mm x 800 mm (open weave)

Contact **Xiros™** on enquiries@xiros.co.uk for further details or visit our website at www.xiros.co.uk

Developed and manufactured by

Xiros™ Ltd

Springfield House, Whitehouse Lane, Leeds LS19 7UE

Tel. +44 (0) 113 238 7202 Fax. +44 (0) 113 238 7201

enquiries@xiros.co.uk

www.xiros.co.uk

Xiros Limited, registered in England No. 1664824.

Xiros is a trademark of Xiros. All rights reserved. © Xiros™ 2024. Worldwide patents and patents pending.

Get involved

At Xiros we are constantly striving to improve outcomes for patients, and we want to hear from you. Help us create the future of sports medicine and share your opinions and expertise with us at:

www.xiros.co.uk